



Human MEPCE blocking peptide (CDBP0600)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Bin3 antibody
Antigen Description	MEPCE (methylphosphate capping enzyme) is a protein-coding gene. Diseases associated with MEPCE include pthirus pubis infestation, and chronic myeloid leukemia. GO annotations related to this gene include S-adenosylmethionine-dependent methyltransferase activity and RNA methyltransferase activity.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MEPCE methylphosphate capping enzyme [Homo sapiens]
Official Symbol	MEPCE
Synonyms	MEPCE; methylphosphate capping enzyme; BCDIN3, bin3, bicoid interacting 3, homolog (Drosophila); 7SK snRNA methylphosphate capping enzyme; FLJ20257; MePCE; bin3 homolog; bin3, bicoid-interacting 3, homolog; bicoid-interacting protein 3 homolog; BCDIN3;
Entrez Gene ID	56257

mRNA Refseq	NM_001194990
Protein Refseq	NP_001181919
UniProt ID	Q7L2J0
Chromosome Location	7q22.1
Function	methyltransferase activity; transferase activity;