



Human PTRH2 blocking peptide (CDBP0603)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-Bit1 antibody
Antigen Description	PTRH2 (peptidyl-tRNA hydrolase 2) is a protein-coding gene. Diseases associated with PTRH2 include esophagitis, and among its related super-pathways are Apoptosis. GO annotations related to this gene include protein binding and aminoacyl-tRNA hydrolase activity.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	PTRH2 peptidyl-tRNA hydrolase 2 [Homo sapiens]
Official Symbol	PTRH2
Synonyms	PTRH2; peptidyl-tRNA hydrolase 2; peptidyl-tRNA hydrolase 2, mitochondrial; Bcl 2 inhibitor of transcription; BIT1; CGI 147; PTH2; PTH 2; bcl-2 inhibitor of transcription 1; CGI-147;

FLJ32471;

Entrez Gene ID	51651
mRNA Refseq	NM_016077
Protein Refseq	NP_057161
UniProt ID	Q9Y3E5
Chromosome Location	17q23.2
Function	aminoacyl-tRNA hydrolase activity; hydrolase activity;