



Human BMP3 blocking peptide (CDBP0611)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-BMP3 antibody
Antigen Description	BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	BMP3 bone morphogenetic protein 3 [Homo sapiens]
Official Symbol	BMP3
Synonyms	BMP3; bone morphogenetic protein 3; bone morphogenetic protein 3 (osteogenic); osteogenin; BMP-3; bone morphogenetic protein-3; bone morphogenetic protein 3A; BMP-3A;
Entrez Gene ID	651
mRNA Refseq	NM_001201

Protein Refseq	NP_001192
UniProt ID	P12645
Chromosome Location	4q21
Pathway	Adipogenesis, organism-specific biosystem;
Function	BMP receptor binding; cytokine activity; growth factor activity; receptor binding;