



Human BORA blocking peptide (CDBP0623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-Bora antibody
Antigen Description	BORA is an activator of the protein kinase Aurora A (AURKA; MIM 603072), which is required for centrosome maturation, spindle assembly, and asymmetric protein localization during mitosis.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	BORA bora, aurora kinase A activator [Homo sapiens]
Official Symbol	BORA
Synonyms	BORA; bora, aurora kinase A activator; C13orf34, chromosome 13 open reading frame 34; protein aurora borealis; FLJ22624; hsBora; aurora borealis; C13orf34; RP11-342J4.2;

Entrez Gene ID	79866
mRNA Refseq	NM_024808
Protein Refseq	NP_079084
UniProt ID	Q6PGQ7
Chromosome Location	13q22.1
Pathway	PLK1 signaling events, organism-specific biosystem;
Function	protein binding; protein kinase binding;