



Human SLC25A47 blocking peptide (CDBP0640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-C14orf68 antibody
Antigen Description	SLC25A47 (solute carrier family 25, member 47) is a protein-coding gene. Diseases associated with SLC25A47 include hepatocellular carcinoma, and hepatoblastoma. An important paralog of this gene is SLC25A45.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SLC25A47 solute carrier family 25, member 47 [Homo sapiens (human)]
Official Symbol	SLC25A47
Synonyms	SLC25A47; solute carrier family 25, member 47; HDMCP; C14orf68; HMFN1655; solute carrier family 25 member 47; HCC-down-regulated mitochondrial carrier protein; hepatocellular carcinoma-downregulated mitochondrial carrier protein; hepatocellular carcinoma down-

regulated mitochondrial carrier protein;

Entrez Gene ID	283600
mRNA Refseq	NM_207117.2
Protein Refseq	NP_997000.2
UniProt ID	Q6Q0C1
Chromosome Location	14q32.2