



Human LRRN4 blocking peptide (CDBP0642)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-C20orf75 antibody
Antigen Description	LRRN4 (leucine rich repeat neuronal 4) is a protein-coding gene. Diseases associated with LRRN4 include neuronitis. An important paralog of this gene is LRRN4CL.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	LRRN4 leucine rich repeat neuronal 4 [Homo sapiens (human)]
Official Symbol	LRRN4
Synonyms	LRRN4; leucine rich repeat neuronal 4; NLRR4; NLRR-4; C20orf75; dJ1056H1.1; leucine-rich repeat neuronal protein 4; neuronal leucine-rich repeat protein 4;
Entrez Gene ID	164312
mRNA Refseq	NM_152611.4

Protein Refseq	NP_689824.2
UniProt ID	Q8WUT4
Chromosome Location	20p12.3