



Human CABP7 blocking peptide (CDBP0647)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-CaBP7 antibody
Antigen Description	CABP7 (calcium binding protein 7) is a protein-coding gene. Diseases associated with CABP7 include retinitis. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CALML3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	CABP7 calcium binding protein 7 [Homo sapiens]
Official Symbol	CABP7
Synonyms	CABP7; calcium binding protein 7; calcium-binding protein 7; MGC57793; calneuron 2; calneuron-2; calneuron II; CALN2;

Entrez Gene ID	164633
mRNA Refseq	NM_182527
Protein Refseq	NP_872333
UniProt ID	Q86V35
Chromosome Location	22q12.2
Function	calcium ion binding;