



Human CCDC69 blocking peptide (CDBP0717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-CCDC69 antibody
Antigen Description	CCDC69 (coiled-coil domain containing 69) is a protein-coding gene. An important paralog of this gene is MTUS1.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	CCDC69 coiled-coil domain containing 69 [Homo sapiens]
Official Symbol	CCDC69
Synonyms	CCDC69; coiled-coil domain containing 69; coiled-coil domain-containing protein 69; DKFZP434C171; FLJ13705; DKFZp434C171;

Entrez Gene ID	26112
mRNA Refseq	NM_015621
Protein Refseq	NP_056436
UniProt ID	A6NI79
Chromosome Location	5q33.1