



Cdc2 blocking peptide (CDBP0739)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Cdc2 peptide
Antigen Description	Transcription of the human cdc2 gene is cell cycle regulated and restricted to proliferating cells. cdc2 transcription is high in S and G2 phases of the cell cycle but low in G1.
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	0.2 mg/ml
Size	100 µg
Buffer	PBS with 100ug BSA 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Keep as concentrated solution, aliquot and store at 4°C.

GENE INFORMATION

Gene Name	cdc2 cell division control protein 2-like protein [Ictalurus punctatus]
Official Symbol	Cdc2
Synonyms	CDC2; cell division control protein 2-like protein;
Entrez Gene ID	100528197
mRNA Refseq	NM_001200467

