



# Human CDC37L1 blocking peptide (CDBP0740)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | Blocking/Immunizing peptide for anti-CDC37L1/HARC antibody                                                                                                                                                          |
| <b>Antigen Description</b> | CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding (Scholz et al., 2001 (PubMed 11413142)). |
| <b>Species</b>             | Human                                                                                                                                                                                                               |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                                                        |
| <b>Applications</b>        | Apuri, BL, ELISA                                                                                                                                                                                                    |
| <b>Format</b>              | Lyophilized powder                                                                                                                                                                                                  |
| <b>Size</b>                | 100 µg                                                                                                                                                                                                              |
| <b>Preservative</b>        | None                                                                                                                                                                                                                |
| <b>Storage</b>             | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                     |

## GENE INFORMATION

|                        |                                                                                                                                                                                                                                                                                                                                              |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">CDC37L1 cell division cycle 37 homolog (S. cerevisiae)-like 1 [ Homo sapiens ]</a>                                                                                                                                                                                                                                               |
| <b>Official Symbol</b> | CDC37L1                                                                                                                                                                                                                                                                                                                                      |
| <b>Synonyms</b>        | CDC37L1; cell division cycle 37 homolog (S. cerevisiae)-like 1; CDC37 cell division cycle 37 homolog (S. cerevisiae) like 1 , cell division cycle 37 homolog (S. cerevisiae) like 1; hsp90 co-chaperone Cdc37-like 1; CDC37B; FLJ20639; HARC; Hsp90-associating relative of Cdc37; CDC37 cell division cycle 37 homolog-like 1; RP11-6J24.5; |
| <b>Entrez Gene ID</b>  | <a href="#">55664</a>                                                                                                                                                                                                                                                                                                                        |

|                     |                           |
|---------------------|---------------------------|
| mRNA Refseq         | <a href="#">NM_017913</a> |
| Protein Refseq      | <a href="#">NP_060383</a> |
| UniProt ID          | Q7L3B6                    |
| Chromosome Location | 9p24.1                    |