



Human CHCHD3 blocking peptide (CDBP0780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CHCHD3 antibody
Antigen Description	CHCHD3 (coiled-coil-helix-coiled-coil-helix domain containing 3) is a protein-coding gene. Diseases associated with CHCHD3 include lung cancer. GO annotations related to this gene include phosphatase binding and RNA polymerase II transcription regulatory region sequence-specific DNA binding transcription factor activity involved in negative regulation of transcription. An important paralog of this gene is CHCHD6.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CHCHD3 coiled-coil-helix-coiled-coil-helix domain containing 3 [Homo sapiens]
Official Symbol	CHCHD3
Synonyms	CHCHD3; coiled-coil-helix-coiled-coil-helix domain containing 3; coiled-coil-helix-coiled-coil-helix domain-containing protein 3, mitochondrial; FLJ20420; MINOS3; mitochondrial inner membrane organizing system 3; PPP1R22; protein phosphatase 1; regulatory subunit 22;

protein phosphatase 1, regulatory subunit 22;

Entrez Gene ID	54927
mRNA Refseq	NM_017812
Protein Refseq	NP_060282
UniProt ID	Q9NX63
Chromosome Location	7q33
Function	phosphatase binding; protein binding; protein complex scaffold;