



# Human CHFR blocking peptide (CDBP0785)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Blocking/Immunizing peptide for anti-CHFR antibody
<b>Antigen Description</b>	CHFR (checkpoint with forkhead and ring finger domains, E3 ubiquitin protein ligase) is a protein-coding gene. Diseases associated with CHFR include malignant peripheral nerve sheath tumor, and laryngeal squamous cell carcinoma, and among its related super-pathways are DNA Damage. GO annotations related to this gene include ubiquitin-protein ligase activity and nucleotide binding.
<b>Species</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Apuri, BL, ELISA
<b>Format</b>	Lyophilized powder
<b>Size</b>	100 µg
<b>Preservative</b>	None
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">CHFR checkpoint with forkhead and ring finger domains [ Homo sapiens ]</a>
<b>Official Symbol</b>	CHFR
<b>Synonyms</b>	CHFR; checkpoint with forkhead and ring finger domains;
<b>Entrez Gene ID</b>	<a href="#">56732</a>