



Human UBQLN4 blocking peptide (CDBP0808)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-CIP75 antibody
Antigen Description	UBQLN4 (ubiquilin 4) is a protein-coding gene. Diseases associated with UBQLN4 include mumps. GO annotations related to this gene include protein binding and polyubiquitin binding. An important paralog of this gene is UBQLN3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	UBQLN4 ubiquilin 4 [Homo sapiens]
Official Symbol	UBQLN4
Synonyms	UBQLN4; ubiquilin 4; C1orf6, chromosome 1 open reading frame 6; ubiquilin-4; A1U; ataxin 1 ubiquitin like interacting protein; UBIN; connexin43-interacting protein of 75 kDa; ataxin-1

interacting ubiquitin-like protein; ataxin-1 ubiquitin-like interacting protein; ataxin-1 ubiquitin-like-interacting protein A1U; A1Up; CIP75; C1orf6;

Entrez Gene ID	56893
mRNA Refseq	NM_020131
Protein Refseq	NP_064516
UniProt ID	Q9NRR5
Chromosome Location	1q21
Pathway	Protein processing in endoplasmic reticulum, organism-specific biosystem; Protein processing in endoplasmic reticulum, conserved biosystem;
Function	polyubiquitin binding; protein binding;