



Human CLMP blocking peptide (CDBP0827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CLMP antibody
Antigen Description	The CTX (see VSIG2, MIM 606011) family of proteins, including ASAM, are type I transmembrane proteins within the Ig superfamily that localize to junctional complexes between endothelial and epithelial cells and may play a role in cell-cell adhesion (Raschperger et al., 2004 (PubMed 14573622)).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CLMP CXADR-like membrane protein [Homo sapiens]
Official Symbol	CLMP
Synonyms	CLMP; CXADR-like membrane protein; ACAM; adipocyte adhesion molecule; adipocyte specific adhesion molecule; ASAM; coxsackie and adenovirus receptor like membrane protein; FLJ22415; CAR-like membrane protein; adipocyte-specific adhesion molecule; coxsackie- and adenovirus receptor-like membrane protein;

Entrez Gene ID	79827
mRNA Refseq	NM_024769
Protein Refseq	NP_079045
UniProt ID	Q9H6B4
Chromosome Location	11q24