



Human CLPB blocking peptide (CDBP0829)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CLPB antibody
Antigen Description	CLPB (ClpB caseinolytic peptidase B homolog (E. coli)) is a protein-coding gene. Diseases associated with CLPB include tularemia, and gastroenteritis. GO annotations related to this gene include nucleoside-triphosphatase activity and molecular_function.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	CLPB ClpB caseinolytic peptidase B homolog (E. coli) [Homo sapiens]
Official Symbol	CLPB
Synonyms	CLPB; ClpB caseinolytic peptidase B homolog (E. coli); caseinolytic peptidase B protein homolog; FLJ13152; HSP78; SKD3; suppressor of potassium transport defect 3;
Entrez Gene ID	81570
mRNA Refseq	NM_030813

Protein Refseq	NP_110440
UniProt ID	Q9H078
Chromosome Location	11q13.4
Function	ATP binding; molecular_function; nucleoside-triphosphatase activity; nucleotide binding;