



Human IFT74 blocking peptide (CDBP0833)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-CMG1/CCDC2/IFT74 antibody
Antigen Description	IFT74 (intraflagellar transport 74 homolog (Chlamydomonas)) is a protein-coding gene. Diseases associated with IFT74 include frontotemporal dementia, and dementia. GO annotations related to this gene include chromatin binding.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	IFT74 intraflagellar transport 74 homolog (Chlamydomonas) [Homo sapiens (human)]
Official Symbol	IFT74
Synonyms	IFT74; intraflagellar transport 74 homolog (Chlamydomonas); CMG1; CCDC2; CMG-1; intraflagellar transport protein 74 homolog; coiled-coil domain containing 2; capillary morphogenesis protein 1; capillary morphogenesis gene 1 protein; coiled-coil domain-containing protein 2;
Entrez Gene ID	80173

mRNA Refseq	NM_001099222.1
Protein Refseq	NP_001092692.1
UniProt ID	A0PJM7
Chromosome Location	9p21.2
Function	beta-tubulin binding; chromatin binding;