



Human COG7 blocking peptide (CDBP0841)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-COG7 antibody
Antigen Description	The protein encoded by this gene resides in the golgi, and constitutes one of the 8 subunits of the conserved oligomeric Golgi (COG) complex, which is required for normal golgi morphology and localization. Mutations in this gene are associated with the congenital disorder of glycosylation type IIe
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	COG7 component of oligomeric golgi complex 7 [Homo sapiens]
Official Symbol	COG7
Synonyms	COG7; component of oligomeric golgi complex 7; conserved oligomeric Golgi complex subunit 7; COG complex subunit 7; CDG2E;
Entrez Gene ID	91949

mRNA Refseq	NM_153603
Protein Refseq	NP_705831
UniProt ID	P83436
Chromosome Location	16p12.2
Function	protein binding;