



# Human CSMD3 blocking peptide (CDBP0892)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Blocking/Immunizing peptide for anti-CSMD3 antibody
<b>Antigen Description</b>	CSMD3 (CUB and Sushi multiple domains 3) is a protein-coding gene. Diseases associated with CSMD3 include exostosis, and familial colorectal cancer. An important paralog of this gene is CSMD1.
<b>Species</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Apuri, BL, ELISA
<b>Format</b>	Lyophilized powder
<b>Size</b>	100 µg
<b>Preservative</b>	None
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">CSMD3 CUB and Sushi multiple domains 3 [ Homo sapiens (human) ]</a>
<b>Official Symbol</b>	CSMD3
<b>Synonyms</b>	CSMD3; CUB and Sushi multiple domains 3; CUB and sushi domain-containing protein 3; CUB and sushi multiple domains protein 3;
<b>Entrez Gene ID</b>	<a href="#">114788</a>
<b>mRNA Refseq</b>	<a href="#">NM_052900.2</a>

<b>Protein Refseq</b>	<a href="#">NP_443132.3</a>
<b>UniProt ID</b>	Q7Z407
<b>Chromosome Location</b>	8q23.3