



Human DBF4B blocking peptide (CDBP0967)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DBF4B antibody
Antigen Description	This gene encodes a regulator of the cell division cycle 7 homolog (<i>S. cerevisiae</i>) protein, a serine-threonine kinase which links cell cycle regulation to genome duplication. This protein localizes to the nucleus and, in complex with the cell division cycle 7 homolog (<i>S. cerevisiae</i>) protein, may facilitate M phase progression. Alternative splicing results in multiple transcript variants.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DBF4B DBF4 homolog B (<i>S. cerevisiae</i>) [<i>Homo sapiens</i>]
Official Symbol	DBF4B
Synonyms	DBF4B; DBF4 homolog B (<i>S. cerevisiae</i>); protein DBF4 homolog B; ASKL1; chifb; chiffon homolog B (<i>Drosophila</i>); DRF1; FLJ13087; ZDBF1B; zinc finger; DBF type containing 1B; chiffon homolog B; ASK-like protein 1; Dbf4-related factor 1; zinc finger, DBF-type containing

1B; activator of S-phase kinase-like protein 1; CHIFB; MGC15009;

Entrez Gene ID	80174
mRNA Refseq	NM_025104
Protein Refseq	NP_079380
UniProt ID	Q8NFT6
Chromosome Location	17q21.31
Function	metal ion binding; nucleic acid binding; zinc ion binding;