



Human DCUN1D1 blocking peptide (CDBP0980)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DCUN1D1 antibody
Antigen Description	DCUN1D1 (DCN1, defective in cullin neddylation 1, domain containing 1) is a protein-coding gene. Diseases associated with DCUN1D1 include squamous cell carcinoma, and adrenocortical carcinoma. GO annotations related to this gene include protein binding. An important paralog of this gene is DCUN1D4.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DCUN1D1 DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae) [Homo sapiens]
Official Symbol	DCUN1D1
Synonyms	DCUN1D1; DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae); DCN1-

like protein 1; DCUN1L1; RP42; SCCRO; SCRO; Tes3; RP42 homolog; DCUN1 domain-containing protein 1; squamous cell carcinoma-related oncogene; defective in cullin neddylation protein 1-like protein 1; DCNL1;

Entrez Gene ID	54165
mRNA Refseq	NM_020640
Protein Refseq	NP_065691
UniProt ID	Q96GG9
Chromosome Location	3q26.3
Function	protein binding;