



# Human DDA1 blocking peptide (CDBP0982)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| Product Overview    | Blocking/Immunizing peptide for anti-DDA1 (aa32-42) antibody |
| Antigen Description | DDA1 (DET1 and DDB1 associated 1) is a protein-coding gene.  |
| Species             | Human                                                        |
| Conjugate           | Unconjugated                                                 |
| Applications        | Apuri, BL, ELISA                                             |
| Format              | Lyophilized powder                                           |
| Size                | 100 µg                                                       |
| Preservative        | None                                                         |
| Storage             | Shipped at ambient temperature, store at -20°C.              |

## GENE INFORMATION

|                 |                                                                                                                                                                                                                                        |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">DDA1 DET1 and DDB1 associated 1 [ Homo sapiens ]</a>                                                                                                                                                                       |
| Official Symbol | DDA1                                                                                                                                                                                                                                   |
| Synonyms        | DDA1; DET1 and DDB1 associated 1; C19orf58, chromosome 19 open reading frame 58; DET1- and DDB1-associated protein 1; MGC2594; PCIA1; PCIA-1; cross-immune reaction antigen PCIA1; placenta cross-immune reaction antigen 1; C19orf58; |
| Entrez Gene ID  | <a href="#">79016</a>                                                                                                                                                                                                                  |
| mRNA Refseq     | <a href="#">NM_024050</a>                                                                                                                                                                                                              |

|                            |                           |
|----------------------------|---------------------------|
| <b>Protein Refseq</b>      | <a href="#">NP_076955</a> |
| <b>UniProt ID</b>          | Q9BW61                    |
| <b>Chromosome Location</b> | 19p13.11                  |