



Human DGAT1 blocking peptide (CDBP0999)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DGAT1 (Internal) antibody
Antigen Description	This gene encodes an multipass transmembrane protein that functions as a key metabolic enzyme. The encoded protein catalyzes the conversion of diacylglycerol and fatty acyl CoA to triacylglycerol. This enzyme can also transfer acyl CoA to retinol. Activity of this protein may be associated with obesity and other metabolic diseases.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DGAT1 diacylglycerol O-acyltransferase 1 [Homo sapiens]
Official Symbol	DGAT1
Synonyms	DGAT1; diacylglycerol O-acyltransferase 1; diacylglycerol O acyltransferase homolog 1 (mouse); ARGPI; DGAT; ACAT related gene product 1; ACAT-related gene product 1; diglyceride acyltransferase; diacylglycerol O-acyltransferase homolog 1; acyl coenzyme A:cholesterol acyltransferase related gene 1;

Entrez Gene ID	8694
mRNA Refseq	NM_012079
Protein Refseq	NP_036211
UniProt ID	O75907
Chromosome Location	8q24.3
Pathway	Fat digestion and absorption, organism-specific biosystem; Fat digestion and absorption, conserved biosystem; Fatty acid, triacylglycerol, and ketone body metabolism, organism-specific biosystem; Glycerolipid metabolism, organism-specific biosystem; Glycerolipid metabolism, conserved biosystem; Metabolic pathways, organism-specific biosystem; Metabolism, organism-specific biosystem;
Function	2-acylglycerol O-acyltransferase activity; O-acyltransferase activity; diacylglycerol O-acyltransferase activity; diacylglycerol binding; fatty acid binding; transferase activity, transferring acyl groups;