



Human DIXDC1 blocking peptide (CDBP1011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DIXDC1 antibody
Antigen Description	The protein encoded by this gene is a positive regulator of the Wnt signaling pathway. The encoded protein is found associated with gamma tubulin at the centrosome. Alternative splicing results in multiple transcript variants encoding different isoforms.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DIXDC1 DIX domain containing 1 [Homo sapiens]
Official Symbol	DIXDC1
Synonyms	DIXDC1; DIX domain containing 1; dixin; KIAA1735; coiled-coil-DIX1; coiled-coil protein DIX1; DIX domain-containing protein 1; CCD1;
Entrez Gene ID	85458
mRNA Refseq	NM_001037954

Protein Refseq	NP_001033043
UniProt ID	Q155Q3
Chromosome Location	11q23.1
Function	actin binding; gamma-tubulin binding; mitogen-activated protein kinase kinase kinase binding; protein binding; protein domain specific binding; signal transducer activity;