



Human DTX3 blocking peptide (CDBP1063)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-DTX3 antibody
Antigen Description	DTX3 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 (PubMed 12670957)).(supplied by OMIM, Nov 2009)
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	DTX3 deltex homolog 3 (Drosophila) [Homo sapiens]
Official Symbol	DTX3
Synonyms	DTX3; deltex homolog 3 (Drosophila); deltex 3 homolog (Drosophila); protein deltex-3; FLJ34766; RNF154; deltex3; deltex 3 homolog; RING finger protein 154; MGC138863; MGC138864;
Entrez Gene ID	196403
mRNA Refseq	NM_178502

Protein Refseq	NP_848597
UniProt ID	Q8N9I9
Chromosome Location	12q13.2
Pathway	Notch signaling pathway, organism-specific biosystem; Notch signaling pathway, conserved biosystem;
Function	metal ion binding; zinc ion binding;