



Human ELMO3 blocking peptide (CDBP1114)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ELMO3 antibody
Antigen Description	The protein encoded by this gene is similar to a <i>C. elegans</i> protein that functions in phagocytosis of apoptotic cells and in cell migration. Other members of this small family of engulfment and cell motility (ELMO) proteins have been shown to interact with the dedicator of cyto-kinesis 1 protein to promote phagocytosis and effect cell shape changes.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ELMO3 engulfment and cell motility 3 [Homo sapiens]
Official Symbol	ELMO3
Synonyms	ELMO3; engulfment and cell motility 3; engulfment and cell motility 3 (ced 12 homolog, <i>C. elegans</i>); engulfment and cell motility protein 3; CED 12; CED12; ELMO 3; FLJ13824; ced-12 homolog 3; CED-12; ELMO-3;
Entrez Gene ID	79767

mRNA Refseq	NM_024712
Protein Refseq	NP_078988
UniProt ID	Q96BJ8
Chromosome Location	16q22.1
Function	SH3 domain binding; binding;