



Human EML2 blocking peptide (CDBP1119)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-EML2/ELP70 antibody
Antigen Description	EML2 (echinoderm microtubule associated protein like 2) is a protein-coding gene. Diseases associated with EML2 include usher syndrome. GO annotations related to this gene include catalytic activity and receptor binding. An important paralog of this gene is EML3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	EML2 echinoderm microtubule associated protein like 2 [Homo sapiens]
Official Symbol	EML2
Synonyms	EML2; echinoderm microtubule associated protein like 2; echinoderm microtubule-associated protein-like 2; echinoderm MT associated protein (EMAP) like protein 70; ELP70; EMAP 2; EMAP2; microtubule associated protein like echinoderm EMAP; microtubule-associated protein like echinoderm EMAP; echinoderm MT-associated protein (EMAP)-like protein 70; EMAP-2;
Entrez Gene ID	24139

mRNA Refseq	NM_001193268
Protein Refseq	NP_001180197
UniProt ID	O95834
Chromosome Location	19q13.32
Function	catalytic activity;