



Human EML5 blocking peptide (CDBP1120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-EML5 antibody
Antigen Description	EML5 (echinoderm microtubule associated protein like 5) is a protein-coding gene. Diseases associated with EML5 include usher syndrome. An important paralog of this gene is EML3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	EML5 echinoderm microtubule associated protein like 5 [Homo sapiens (human)]
Official Symbol	EML5
Synonyms	EML5; echinoderm microtubule associated protein like 5; EMAP-2; echinoderm microtubule-associated protein-like 5; EMAP-5;
Entrez Gene ID	161436
mRNA Refseq	NM_183387.2

Protein Refseq	NP_899243.1
UniProt ID	Q05BV3
Chromosome Location	14q31.3
Function	catalytic activity; molecular_function;