



Human EMP3 blocking peptide (CDBP1121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-EMP3 antibody
Antigen Description	EMP3 (epithelial membrane protein 3) is a protein-coding gene. Diseases associated with EMP3 include pheochromocytoma, and keratoconus. An important paralog of this gene is EMP1.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	EMP3 epithelial membrane protein 3 [Homo sapiens]
Official Symbol	EMP3
Synonyms	EMP3; epithelial membrane protein 3; YMP; EMP-3; HNMP-1; hematopoietic neural membrane protein 1;
Entrez Gene ID	2014
mRNA Refseq	NM_001425

Protein Refseq	NP_001416
UniProt ID	P54852
Chromosome Location	19q13.3