



Human HECW2 blocking peptide (CDBP1125)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-EndoNEDD4/NEDL2 antibody
Antigen Description	HECW2 (HECT, C2 and WW domain containing E3 ubiquitin protein ligase 2) is a protein-coding gene. Diseases associated with HECW2 include liddle syndrome, and endotheliitis. GO annotations related to this gene include ubiquitin-protein ligase activity. An important paralog of this gene is NEDD4L.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HECW2 HECT, C2 and WW domain containing E3 ubiquitin protein ligase 2 [Homo sapiens (human)]
Official Symbol	HECW2
Synonyms	HECW2; HECT, C2 and WW domain containing E3 ubiquitin protein ligase 2; NEDL2; E3 ubiquitin-protein ligase HECW2; NEDD4-related E3 ubiquitin ligase NEDL2; NEDD4-like E3 ubiquitin-protein ligase 2; HECT, C2 and WW domain-containing protein 2;

Entrez Gene ID	57520
mRNA Refseq	NM_020760.1
Protein Refseq	NP_065811.1
UniProt ID	Q9P2P5
Chromosome Location	2q32.3
Function	ubiquitin-protein ligase activity;