



Human EPHB4 blocking peptide (CDBP1140)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-EPHB4 (aa 566 - 578) antibody
Antigen Description	Ephrin receptors and their ligands, the ephrins, mediate numerous developmental processes, particularly in the nervous system. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. The Eph family of receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. Ephrin receptors make up the largest subgroup of the receptor tyrosine kinase (RTK) family. The protein encoded by this gene binds to ephrin-B2 and plays an essential role in vascular development. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	EPHB4 EPH receptor B4 [Homo sapiens]
Official Symbol	EPHB4

Synonyms	EPHB4; EPH receptor B4; EphB4 , HTK; ephrin type-B receptor 4; Tyro11; ephrin receptor EphB4; soluble EPHB4 variant 1; soluble EPHB4 variant 2; soluble EPHB4 variant 3; hepatoma transmembrane kinase; tyrosine-protein kinase TYRO11; tyrosine-protein kinase receptor HTK; HTK; MYK1; TYRO11;
Entrez Gene ID	2050
mRNA Refseq	NM_004444
Protein Refseq	NP_004435
UniProt ID	P54760
Chromosome Location	7q22
Pathway	Axon guidance, organism-specific biosystem; Axon guidance, conserved biosystem; EPHB forward signaling, organism-specific biosystem; Ephrin B reverse signaling, organism-specific biosystem; EphrinB-EPHB pathway, organism-specific biosystem;
Function	ATP binding; ephrin receptor activity; ephrin receptor activity; nucleotide binding; protein binding; protein tyrosine kinase activity; receptor activity; transmembrane receptor protein tyrosine kinase activity;