



Human ETV2 blocking peptide (CDBP1165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ETV2 antibody
Antigen Description	ETV2 (ets variant 2) is a protein-coding gene. Diseases associated with ETV2 include albright's hereditary osteodystrophy, and tonsillitis. GO annotations related to this gene include sequence-specific DNA binding and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is ELK4.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ETV2 ets variant 2 [Homo sapiens]
Official Symbol	ETV2
Synonyms	ETV2; ets variant 2; ets variant gene 2; ETS translocation variant 2; ER71; ets-related protein 71; ETSRP71; MGC129834; MGC129835;
Entrez Gene ID	2116

mRNA Refseq	NM_014209
Protein Refseq	NP_055024
UniProt ID	O00321
Chromosome Location	19q13.11
Function	sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;