



Human FBXL12 blocking peptide (CDBP1203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FBXL12 antibody
Antigen Description	Members of the F-box protein family, such as FBXL12, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FBXL12 F-box and leucine-rich repeat protein 12 [Homo sapiens]
Official Symbol	FBXL12
Synonyms	FBXL12; F-box and leucine-rich repeat protein 12; F-box/LRR-repeat protein 12; Fbl12; FLJ20188; F-box protein FBL12;
Entrez Gene ID	54850

mRNA Refseq	NM_017703
Protein Refseq	NP_060173
UniProt ID	Q9NXK8
Chromosome Location	19p13.2
Function	protein binding;