



Human FBXO6 blocking peptide (CDBP1211)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FBXO6 antibody
Antigen Description	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class, and its C-terminal region is highly similar to that of rat NFB42 (neural F Box 42 kDa) which may be involved in the control of the cell cycle.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FBXO6 F-box protein 6 [Homo sapiens]
Official Symbol	FBXO6

Synonyms	FBXO6; F-box protein 6; F box only protein 6; F-box only protein 6; FBG2; FBS2; FBX6; Fbx6b; F-box protein FBG2; F-box protein Fbx6; F-box/G-domain protein 2; F-box protein that recognizes sugar chains 2;
Entrez Gene ID	26270
mRNA Refseq	NM_018438
Protein Refseq	NP_060908
UniProt ID	Q9NRD1
Chromosome Location	1p36.22
Pathway	Adaptive Immune System, organism-specific biosystem; Antigen processing: Ubiquitination & Proteasome degradation, organism-specific biosystem; Association of TriC/CCT with target proteins during biosynthesis, organism-specific biosystem; Chaperonin-mediated protein folding, organism-specific biosystem; Class I MHC mediated antigen processing & presentation, organism-specific biosystem;
Function	carbohydrate binding; glycoprotein binding; protein binding; ubiquitin-protein ligase activity;