



Human FEM1A blocking peptide (CDBP1214)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FEM1A/APRAP antibody
Antigen Description	FEM1A (fem-1 homolog a (C. elegans)) is a protein-coding gene. Diseases associated with FEM1A include polycystic ovary syndrome, and rhabdomyosarcoma. GO annotations related to this gene include ubiquitin-protein ligase activity and EP4 subtype prostaglandin E2 receptor binding. An important paralog of this gene is TANC2.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FEM1A fem-1 homolog a (C. elegans) [Homo sapiens]
Official Symbol	FEM1A
Synonyms	FEM1A; fem-1 homolog a (C. elegans); protein fem-1 homolog A; FEM1-alpha; prostaglandin E receptor 4-associated protein; EPRAP; DKFZp762M136;
Entrez Gene ID	55527

mRNA Refseq	NM_018708
Protein Refseq	NP_061178
UniProt ID	Q9BSK4
Chromosome Location	19p13.3
Function	EP4 subtype prostaglandin E2 receptor binding; binding; ubiquitin-protein ligase activity;