



Human FEZ2 blocking peptide (CDBP1218)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-FEZ2 antibody
Antigen Description	This gene is an ortholog of the C. elegans unc-76 gene, which is necessary for normal axonal bundling and elongation within axon bundles. Other orthologs include the rat gene that encodes zygin II, which can bind to synaptotagmin. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	FEZ2 fasciculation and elongation protein zeta 2 (zygin II) [Homo sapiens (human)]
Official Symbol	FEZ2
Synonyms	FEZ2; fasciculation and elongation protein zeta 2 (zygin II); HUM3CL; fasciculation and elongation protein zeta-2; zygin 2; zygin-2; zygin II; pre-T/NK cell associated protein (3CI);

Entrez Gene ID	9637
mRNA Refseq	NM_001042548.1
Protein Refseq	NP_001036013.1
UniProt ID	Q9UHY8
Chromosome Location	2p21
Function	protein binding;