



Human FMR1 blocking peptide (CDBP1248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FMR1 antibody
Antigen Description	The protein encoded by this gene binds RNA and is associated with polysomes. The encoded protein may be involved in mRNA trafficking from the nucleus to the cytoplasm. A trinucleotide repeat (CGG) in the 5' UTR is normally found at 6-53 copies, but an expansion to 55-230 repeats is the cause of fragile X syndrome. Expansion of the trinucleotide repeat may also cause one form of premature ovarian failure (POF1). Multiple alternatively spliced transcript variants that encode different protein isoforms and which are located in different cellular locations have been described for this gene. [provided by RefSeq, May 2010]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FMR1 fragile X mental retardation 1 [Homo sapiens (human)]
Official Symbol	FMR1
Synonyms	FMR1; fragile X mental retardation 1; POF; FMRP; POF1; FRAXA; fragile X mental retardation

protein 1;

Entrez Gene ID	2332
mRNA Refseq	NM_001185075.1
Protein Refseq	NP_001172004.1
UniProt ID	G8JLE9
Chromosome Location	Xq27.3
Pathway	RNA transport, organism-specific biosystem; RNA transport, conserved biosystem;
Function	RNA binding; mRNA binding; poly(A) RNA binding; protein binding;