



Human FNIP2 blocking peptide (CDBP1250)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-FNIP2 antibody
Antigen Description	FNIP2 (folliculin interacting protein 2) is a protein-coding gene. Diseases associated with FNIP2 include oxyphilic adenoma, and pneumothorax. GO annotations related to this gene include protein binding. An important paralog of this gene is FNIP1.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	FNIP2 folliculin interacting protein 2 [Homo sapiens]
Official Symbol	FNIP2
Synonyms	FNIP2; folliculin interacting protein 2; folliculin-interacting protein 2; FNIP1; KIAA1450; FNIP1-like protein; O6-methylguanine-induced apoptosis 1 protein; MAPO1;

Entrez Gene ID	57600
mRNA Refseq	NM_020840
Protein Refseq	NP_065891
UniProt ID	Q9P278
Chromosome Location	4q32.1
Function	protein binding;