



Mouse Foxi2 blocking peptide (CDBP1266)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Foxi2 (mouse) antibody
Antigen Description	FOXI2 (forkhead box I2) is a protein-coding gene. GO annotations related to this gene include double-stranded DNA binding and sequence-specific DNA binding. An important paralog of this gene is FOXF1.
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Foxi2 forkhead box I2 [Mus musculus]
Official Symbol	Foxi2
Synonyms	FOXI2; forkhead box I2; forkhead box protein I2; B130055A05Rik;
Entrez Gene ID	270004
mRNA Refseq	NM_183193

Protein Refseq[NP_899016](#)**Function**

DNA binding; DNA binding, bending; double-stranded DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; sequence-specific distal enhancer binding RNA polymerase II transcription factor activity; transcri