



Human FOXL3 blocking peptide (CDBP1267)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FOXL3/Forkhead box I3 antibody
Antigen Description	FOXL3 (forkhead box I3) is a protein-coding gene. GO annotations related to this gene include double-stranded DNA binding and sequence-specific DNA binding.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FOXL3 forkhead box I3 [Homo sapiens]
Official Symbol	FOXL3
Synonyms	FOXL3; forkhead box I3; forkhead box protein I3;
Entrez Gene ID	344167
mRNA Refseq	NM_001135649
Protein Refseq	NP_001129121

UniProt ID	A8MTJ6
Chromosome Location	2p11.2
Function	DNA binding, bending; double-stranded DNA binding; sequence-specific DNA binding; sequence-specific distal enhancer binding RNA polymerase II transcription factor activity; transcription factor binding;