



Human FSTL1 blocking peptide (CDBP1300)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-FSTL1 antibody
Antigen Description	This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	FSTL1 follistatin-like 1 [Homo sapiens]
Official Symbol	FSTL1
Synonyms	FSTL1; follistatin-like 1; follistatin-related protein 1; FRP; FSL1; follistatin-like protein 1; FLJ50214; FLJ52277;
Entrez Gene ID	11167
mRNA Refseq	NM_007085

Protein Refseq	NP_009016
UniProt ID	Q12841
Chromosome Location	3q13.33
Pathway	Signal Transduction, organism-specific biosystem; Signaling by BMP, organism-specific biosystem;
Function	calcium ion binding; heparin binding;