



Human SYF2 blocking peptide (CDBP1345)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GCIP-interacting protein p29 antibody
Antigen Description	This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SYF2 SYF2 homolog, RNA splicing factor (S. cerevisiae) [Homo sapiens]
Official Symbol	SYF2
Synonyms	SYF2; SYF2 homolog, RNA splicing factor (S. cerevisiae); CBPIN, CCNDBP1 interactor; pre-mRNA-splicing factor SYF2; DKFZp564O2082; fSAP29; functional spliceosome associated protein 29; NTC31; p29; CCNDBP1 interactor; CCNDBP1-interactor; GCIP-interacting protein p29; functional spliceosome-associated protein 29; P29; CBPIN;
Entrez Gene ID	25949

mRNA Refseq	NM_015484
Protein Refseq	NP_056299
UniProt ID	O95926
Chromosome Location	1p36.11
Pathway	Spliceosome, organism-specific biosystem; Spliceosome, conserved biosystem; Spliceosome, 35S U5-snRNP, organism-specific biosystem;