



Human ARHGEF5 blocking peptide (CDBP1354)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GEF5 antibody
Antigen Description	Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ARHGEF5 Rho guanine nucleotide exchange factor (GEF) 5 [Homo sapiens]
Official Symbol	ARHGEF5
Synonyms	ARHGEF5; Rho guanine nucleotide exchange factor (GEF) 5;
Entrez Gene ID	54385