



Human GIPC3 blocking peptide (CDBP1363)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GIPC3 antibody
Antigen Description	The protein encoded by this gene belongs to the GIPC family. Studies in mice suggest that this gene is required for postnatal maturation of the hair bundle and long-term survival of hair cells and spiral ganglion in the ear. Mutations in this gene are associated with autosomal recessive deafness.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GIPC3 GIPC PDZ domain containing family, member 3 [Homo sapiens]
Official Symbol	GIPC3
Synonyms	GIPC3; GIPC PDZ domain containing family, member 3; C19orf64, chromosome 19 open reading frame 64 , deafness, autosomal recessive 15 , deafness, autosomal recessive 72 , DFNB15, DFNB72; PDZ domain-containing protein GIPC3; DFNB95; DFNB15; DFNB72; C19orf64; FLJ40925; DKFZp686J1198;

Entrez Gene ID	126326
mRNA Refseq	NM_133261
Protein Refseq	NP_573568
UniProt ID	Q8TF64
Chromosome Location	19p13.3