



Human GLIS1 blocking peptide (CDBP1373)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GLIS1 antibody
Antigen Description	GLIS1 is a GLI (MIM 165220)-related Kruppel-like zinc finger protein that functions as an activator and repressor of transcription (Kim et al., 2002 [PubMed 12042312]).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GLIS1 GLIS family zinc finger 1 [Homo sapiens]
Official Symbol	GLIS1
Synonyms	GLIS1; GLIS family zinc finger 1; zinc finger protein GLIS1; FLJ36155; GLI-similar 1;
Entrez Gene ID	148979
mRNA Refseq	NM_147193
Protein Refseq	NP_671726

UniProt ID	Q8NBF1
Chromosome Location	1p32.3
Function	DNA binding; metal ion binding; zinc ion binding;