



Human GPR15 blocking peptide (CDBP1408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	GPR15 (N - term) peptide (murine)
Antigen Description	This gene encodes a G protein-coupled receptor that acts as a chemokine receptor for human immunodeficiency virus type 1 and 2. The encoded protein localizes to the cell membrane. [provided by RefSeq, Nov 2012]
Species	Human
Conjugate	Unconjugated
Applications	BL
Concentration	0.2 mg/ml
Size	50 µg
Buffer	Preservative: 0.02% Sodium Azide; Constituents: 0.1% BSA, PBS. pH 7.2
Preservative	0.02% Sodium Azide

GENE INFORMATION

Gene Name	GPR15 G protein-coupled receptor 15 [Homo sapiens (human)]
Official Symbol	GPR15
Synonyms	GPR15; G protein-coupled receptor 15; BOB; G-protein coupled receptor 15; brother of Bonzo;
Entrez Gene ID	2838
mRNA Refseq	NM_005290.2

Protein Refseq	NP_005281.1
UniProt ID	B2R8H6
Chromosome Location	3q11.2-q13.1
Pathway	GPCRs, Class A Rhodopsin-like, organism-specific biosystem;
Function	G-protein coupled receptor activity;