



Human GRIPAP1 blocking peptide (CDBP1426)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-GRASP1/GRIPAP1 antibody
Antigen Description	This gene encodes a guanine nucleotide exchange factor for the Ras family of small G proteins (RasGEF). The encoded protein interacts in a complex with glutamate receptor interacting protein 1 (GRIP1) and plays a role in the regulation of AMPA receptor function.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GRIPAP1 GRIP1 associated protein 1 [Homo sapiens]
Official Symbol	GRIPAP1
Synonyms	GRIPAP1; GRIP1 associated protein 1; GRIP1-associated protein 1; DKFZp434P0630; GRASP 1; GRASP1; KIAA1167; MPMGp800B12492Q3; GRASP-1; MGC126593; MGC126595;
Entrez Gene ID	56850
mRNA Refseq	NM_020137

Protein Refseq	NP_064522
UniProt ID	Q4V328
Chromosome Location	Xp11.23
Function	molecular_function;