



Human HEXIM1 blocking peptide (CDBP1469)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-HEXIM1 antibody
Antigen Description	Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HEXIM1 hexamethylene bis-acetamide inducible 1 [Homo sapiens]
Official Symbol	HEXIM1
Synonyms	HEXIM1; hexamethylene bis-acetamide inducible 1; protein HEXIM1; CLP 1; EDG1; HIS1; MAQ1; HMBA-inducible; menage a quatre 1; cardiac lineage protein 1; menage a quatre protein 1; estrogen down-regulated gene 1 protein; hexamethylene bis-acetamide inducible 1; hexamethylene bisacetamide-inducible protein; hexamethylene bis-acetamide-inducible protein 1; hexamethylene-bis-acetamide-inducible transcript 1; CLP1; FLJ13562;
Entrez Gene ID	10614

mRNA Refseq	NM_006460
Protein Refseq	NP_006451
UniProt ID	O94992
Chromosome Location	17q21.31
Function	cyclin-dependent protein kinase inhibitor activity; protein binding; snRNA binding;