



Human HRNR blocking peptide (CDBP1495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Hornerin antibody
Antigen Description	HRNR (hornerin) is a protein-coding gene. Diseases associated with HRNR include irritant dermatitis, and dermatitis. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is RPTN.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HRNR hornerin [Homo sapiens (human)]
Official Symbol	HRNR
Synonyms	HRNR; hornerin; FLG3; S100A16; S100a18; filaggrin family member 3; intermediate filament-associated protein;
Entrez Gene ID	388697
mRNA Refseq	NM_001009931.2

Protein Refseq	NP_001009931.1
UniProt ID	Q5DT20
Chromosome Location	1q21.3
Function	calcium ion binding;