



Human HPS6 blocking peptide (CDBP1509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-HPS6 antibody
Antigen Description	This intronless gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. This protein interacts with Hermansky-Pudlak syndrome 5 protein. Mutations in this gene are associated with Hermansky-Pudlak syndrome type 6.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HPS6 Hermansky-Pudlak syndrome 6 [Homo sapiens]
Official Symbol	HPS6
Synonyms	HPS6; Hermansky-Pudlak syndrome 6; Hermansky-Pudlak syndrome 6 protein; FLJ22501; ruby-eye protein homolog; Hermansky-Pudlak syndrome-6 protein (HPS6); MGC20522; RP11-302K17.1;
Entrez Gene ID	79803

mRNA Refseq	NM_024747
Protein Refseq	NP_079023
UniProt ID	Q86YV9
Chromosome Location	10q24.32