



Human HUWE1 blocking peptide (CDBP1528)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-HUWE1/LASU1 antibody
Antigen Description	This gene encodes a protein containing a C-terminal HECT (E6AP type E3 ubiquitin protein ligase) domain that functions as an E3 ubiquitin ligase. The encoded protein is required for the ubiquitination and subsequent degradation of the anti-apoptotic protein Mcl1 (myeloid cell leukemia sequence 1 (BCL2-related)). This protein also ubiquitinates the p53 tumor suppressor, core histones, and DNA polymerase beta. Mutations in this gene are associated with Turner type X-linked syndromic mental retardation. [provided by RefSeq, Aug 2013]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HUWE1 HECT, UBA and WWE domain containing 1, E3 ubiquitin protein ligase [Homo sapiens]
Official Symbol	HUWE1
Synonyms	HUWE1; HECT, UBA and WWE domain containing 1, E3 ubiquitin protein ligase; HECT, UBA

and WWE domain containing 1; E3 ubiquitin-protein ligase HUWE1; Ib772; KIAA0312; UREB1; ARF-binding protein 1; URE-binding protein 1; BJ-HCC-24 tumor antigen; large structure of UREB1; HECT domain protein LASU1; Mcl-1 ubiquitin ligase E3; upstream regulatory element-binding protein 1; homologous to E6AP carboxyl terminus homologous protein 9; MULE; LASU1; HECTH9; URE-B1; ARF-BP1; HSPC272;

Entrez Gene ID	10075
mRNA Refseq	NM_031407
Protein Refseq	NP_113584
UniProt ID	Q7Z6Z7
Chromosome Location	Xp11.22
Pathway	Adaptive Immune System, organism-specific biosystem; Antigen processing: Ubiquitination & Proteasome degradation, organism-specific biosystem; Class I MHC mediated antigen processing & presentation, organism-specific biosystem; Immune System, organism-specific biosystem; Ubiquitin mediated proteolysis, organism-specific biosystem;
Function	DNA binding; acid-amino acid ligase activity; binding; ligase activity; protein binding; ubiquitin-protein ligase activity;